



## HILL-ROM ADVANCE

**Bed Model: P1104/P1105/P1115/P1135/  
P1145/P1155/P1165**

Height Overall – Foot End – High Position.....	41 ½" (105.4 cm)
Height Overall – Foot End – Low Position.....	27 ½" (69.8 cm)
Height Overall – Head End – High Position.....	43" (109.2 cm)
Height Overall – Head End – Low Position.....	29" (73.7 cm)
Height – Floor to Top of Surface – High Position.....	31" (78.7 cm)
Height – Floor to Top of Surface – Low Position.....	17" (43.2 cm)
Overall Length (with cane or roller bumpers) – Fully Extended.....	93 ½" (237.5 cm)
Overall Length (with cane or roller bumpers) – Fully Retracted.....	81" (205.7 cm)
Overall Length (without bumpers) – Fully Extended.....	91" (231.1 cm)
Overall Length (without bumpers) – Fully Retracted.....	78 ½" (199.4 cm)
Overall Width (siderails up).....	41 ½" (105.4 cm)
Overall Width (sideguard stored with #818 bumpers).....	36 ½" (92.7 cm)
Total Hilow Travel.....	14" (35.6 cm)
Sleeping Surface.....	36" x 82" (91.4 cm x 208.3 cm)
Bed Frame Clearance.....	7 ½" (19 cm)
Weight – Bed with Casters (model P1104).....	380 lb (172.4 kg)
Weight – Bed with Casters (model P1105).....	385 lb (174.6 kg)
Weight – Bed with Casters (model P1145/P1165).....	445 lb (201.8 kg)
Maximum Patient Weight.....	500 lb (226.8 kg)



**Piedmont Medical Incorporated**

500 W. Kapp St.

Dobson, NC 27017

[www.piedmontmedicalinc.com](http://www.piedmontmedicalinc.com)

Phone: 800-433-3255

Fax: 877-433-3255