

Hill-Rom 8020 TranStar Electric Stretcher

Product Description

- Height Overall – High Position: **35" (88.9 cm)**
- Height Overall – Low Position: **21 ½" (54.6 cm)**
- Overall Length: **83" (210.9 cm)**
- Overall Width (siderails up): **36" (91.4 cm)**
- Upper Frame Width (siderails stored): **30 ¼" (76.8 cm)**
- Overall Siderail Length: **47" (119.4 cm)**
- Siderail Height Above Sleep Surface: **14" (35.6 cm)**
- Mattress Size: **30" x 75" (76.2 cm x 190.5 cm)**
- Caster Size: **8" (20.3 cm)**
- Floor-to-Base Clearance: **3 ½" (8.9 cm)**
- Floor-to-Hydraulic Cylinder Clearance: **2" (5 cm)**
- Base Shroud Width: **28" (71.1 cm)**
- Maximum Back Incline: **90°**
- Maximum Thigh Section Incline: **25°**
- Weight – Stretcher with Casters: **305 – 335 lb. (138.3 – 152 kg)**
- Maximum Patient Weight: **700 lb. (317.5 kg)**

