

# **Automatic Overbed Table**

## **Overbed Table Model: 1010H**

### **TOP:**

Wood grained melamine laminate  
on 3/4" particleboard core. Flush mounted  
T-mold edge. 15"W x 30"L melamine top

### **FRAME:**

H frame chrome plated 16 and 17 gage welded  
tubular steel column and base. 13 gage  
stamped steel support bracket.  
Spring assisted lift (28" to 45"H)

### **CASTERS:**

1 1/2" twin wheel



**Weight Limit: 50 Lbs. evenly distributed over the column.**

**Weight Package: 27 Lbs. Cube 1.1. UPS shippable.**

### **LAMINATE COLOR OPTIONS:**



**Bannister Oak**  
(Light Oak)



**English Oak**  
(Medium Oak)



**Gunstock Walnut**  
(Dark Oak)



**Wild Cherry**