

# Stryker 721 Transport Stretcher

## Product Description

- Overall Length: **83" (210.8 cm)**
- Overall Width (siderails up): **30" (76.2 cm)**
- Height Overall – High Position: **36" (91.4 cm)**
- Height Overall – Low Position: **24 ½" (62.2 cm)**
- Patient Surface: **24" x 76" (61 cm x 193 cm)**
- Siderails: **13" x 55" (33 cm x 137.2 cm)**
- Caster Size: **8" (20.3 cm)**
- Backrest: **0° – 90°**
- Trendelenburg/Reverse Trendelenburg: **±18°**
- Maximum Patient Weight: **500 lb. (226.8 kg)**

