

## Hill-Rom 1115 Advance Bed

### Product Description

- Height Overall – Foot End – High Position **41 ½" (105.4 cm)**
- Height Overall – Foot End – Low Position **27 ½" (69.8 cm)**
- Height Overall – Head End – High Position **43" (109.2 cm)**
- Height Overall – Head End – Low Position **29" (73.7 cm)**
- Height – Floor to Top of Surface – High Position **31" (78.7 cm)**
- Height- Floor to Top of Surface – Low Position **17" (43.2 cm)**
- Overall Length (with cane or roller bumpers) – Fully Extended **93 ½" (237.5 cm)**
- Overall Length (with cane or roller bumpers) – Fully Retracted **81" (205.7 cm)**
- Overall Length (without bumpers) – Fully Extended **91" (231.1 cm)**
- Overall Length (without bumpers) – Fully Retracted **78 ½" (199.4 cm)**
- Overall Width (siderails up) **41 ½" (105.4 cm)**
- Overall Width (sideguard stored with #818 bumpers) **36 ½" (92.7 cm)**
- Total Hi-Low Travel **14" (35.6 cm)**
- Sleeping Surface **36" x 82" (91.4 cm x 208.3 cm)**
- Bed Frame Clearance **7 ½" (19 cm)**
- Weight – Bed with Casters (model P1104) **380 lb. (172.4 kg)**
- Weight – Bed with Casters (model P1105) **385 lb. (174.6 kg)**
- Weight – Bed with Casters (model P1145/P1165) **445 lb. (201.8 kg)**
- Maximum Patient Weight **500 lb. (226.8 kg)**

