

Hill-Rom 1145 Advance Bed

Product Description

- Height Overall – Foot End – High Position: **41 ½" (105.4 cm)**
- Height Overall – Foot End – Low Position: **27 ½" (69.8 cm)**
- Height Overall – Head End – High Position: **43" (109.2 cm)**
- Height Overall – Head End – Low Position: **29" (73.7 cm)**
- Height – Floor to Top of Surface – High Position: **31" (78.7 cm)**
- Height – Floor to Top of Surface – Low Position: **17" (43.2 cm)**
- Overall Length (with cane or roller bumpers) – Fully Extended: **93 ½" (237.5 cm)**
- Overall Length (with cane or roller bumpers) – Fully Retracted: **81" (205.7 cm)**
- Overall Length (without bumpers) – Fully Extended: **91" (231.1 cm)**
- Overall Length (without bumpers) – Fully Retracted: **78 ½" (199.4 cm)**
- Overall Width (siderails up): **41 ½" (105.4 cm)**
- Overall Width (sideguard stored with #818 bumpers): **36 ½" (92.7 cm)**
- Total Hi/low Travel: **14" (35.6 cm)**
- Sleeping Surface: **36" x 82" (91.4 cm x 208.3 cm)**
- DynamicAire Mattress Bed Frame Clearance: **7 ½" (19 cm)**
- Weight – Bed with Casters: **445 lb. (201.8 kg)**
- Maximum Patient Weight: **500 lb. (226.8 kg)**

