

## Hill-Rom 850-852 Centra Bed

### Product Description

- Height Overall – Foot End – High Position: **42" (106.7 cm)**
- Height Overall – Foot End – Low Position: **28" (71.1 cm)**
- Height Overall – Head End – High Position: **43" (109.2 cm)**
- Height Overall – Head End – Low Position: **29" (73.7 cm)**
- Height – Floor to Top of Surface – High Position: **31" (78.7 cm)**
- Height – Floor to Top of Surface – Low Position: **17" (43.2 cm)**
- Overall Length (with cane or roller bumpers) – Fully Extended: **93 ½" (237.5 cm)**
- Overall Length (with cane or roller bumpers) – Fully Retracted: **79 ½" (201.9 cm)**
- Overall Length (without bumpers) – Fully Extended: **91" (231.1 cm)**
- Overall Length (without bumpers) – Fully Retracted: **77" (195.6 cm)**
- Overall Width (siderails up): **41 ½" (105.4 cm)**
- Overall Width (sideguard stored with #818 bumpers): **36 ½" (92.7 cm)**
- Total Hollow Travel: **14" (35.6 cm)**
- Sleeping Surface: **36" x 82" (91.4 cm x 208.3 cm)**
- Weight – Bed with Casters: **346 lb. (156.9 kg)**
- Weight – Bed with Casters and Bed Panels: **373 lb. (169.2 kg)**
- Maximum Patient Weight: **500 lb. (226.8 kg)**
- 4 IV Pole Accommodation
- \*Model 850 / Spring Surface
- \*Model 852 / Hardpan Surface

