

Hill-Rom P1600 Advanta Bed with Scale

Product Description

- Overall Length: **91" (231 cm)**
- Retracted Length: **78 ½" (199.4 cm)**
- Mattress Support Width: **35" (89 cm)**
- Mattress Support Length: **83" (211 cm)**
- Maximum Width (siderails up): **42 ¾" (108.59 cm)**
- Maximum Width (siderails stored): **36" (91 cm)**
- Maximum Headboard Height: **48" (122 cm)**
- Maximum Footboard Height: **45" (114 cm)**
- Maximum Siderail Height: **45" (114 cm)**
- Bed Clearance: **4 ¾" (12.07 cm)**
- Head Section Inclination (maximum): **65°**
- Knee Section Inclination (maximum): **35°**
- Automatic Contour Up (maximum): **15°**
- Automatic Contour Down: **28°**
- Automatic Contour: **10°**
- Hi/low Range: **13 ¼" (33.66 cm)**
- Trendelenburg: **11°**
- Reverse Trendelenburg: **11°**
- Maximum Patient Weight: **500 lb. (226.8 kg)**

