

Hill-Rom P1900 TotalCare with Scale

Product Description

- Height – From Floor to Top of Articulating Deck (bottom of mattress) – Low Position – Between Stepped Sides: **4 ½” – 15 ½” (36.8 – 39.4 cm)**
- Height – From Floor to Top of Articulating Deck (bottom of mattress) – High Position – Between Stepped Sides: **33 ½” – 34 ½” (85.9 – 87.6 cm)**
- Minimum Under Bed Clearance: **1 ¾” (4.45 cm)**
- Minimum Clearance Between Articulating Deck Foot Section and the Floor for any Bed Position, except the Tilt-Table and the Chair Egress Positions with the Foot Section at 90°: **2” (5.1 cm)**
- Overall Width – Siderails Raised: **40” (101.6 cm)**
- Overall Width – Siderail Lowered: **6 ½” (92.7 cm)**
- Overall Width – Distance between Casters: **25 ½” (64.8 cm)**
- Overall Length: **91 ½” (232.4 cm)**
- Overall Length (with roller bumpers): **92 ½” (235 cm)**
- Overall Length (transport shelf installed): **93 ½” (237.5 cm)**
- Overall Foot Section Width: **28 22/25” – 29.13” (73.4 – 74 cm)**
- Fully Extended Length of the Foot Section from Pivot to inside of the Footboard: **26 ¾” – 27 ¼” (67.95 – 69.22 cm)**
- Articulating Deck Length: **81 ¾” – 82 ¼” (207.65 – 208.92 cm)**
- Maximum Head Elevation: **75°**
- Maximum Knee Elevation: **20°**

