

Hill-Rom P3600 Affinity II Birthing Bed

Product Description

- Overall Length: **90" (228.6 cm)**
- Overall Height – Headboard – High Position: **53 ¼" (136 cm)**
- Overall Height – Low Position (without mattress): **18 ½" (47 cm)**
- Overall Height – High Position (without mattress): **5 ¾" (91.44 cm)**
- Total Hollow Travel: **17 ¼" (44.45 cm)**
- Overall Width (siderails down): **34" (86.4 cm)**
- Siderail Length: **28" (71 cm)**
- Clearance Inside Siderails: **38" (96.5 cm)**
- Siderail Height (without mattress): **13" (33 cm)**
- Mattress Size: **34" x 78" (86.4 cm x 198 cm)**
- Caster Size: **6" (15.2 cm)**
- Clearance under Base – High Position: **5" (12.7 cm)**
- Clearance under Base – Low Position: **2 ½" (6.35 cm)**
- Trendelenburg Angle Inclination: **8°**
- Foot Section – Maximum Lift Weight: **400 lb. (181.6 kg)**
- Head Section – Maximum Lift Weight: **200 lb. (91 kg)**
- Maximum Patient Weight: **500 lb. (226.8 kg)**
- ***Model P3600 / Affinity Bed with V-Cut Seat**
- ***Model P3601 / Affinity Bed with U-Cut Seat**

