

Hill-Rom P3700 Affinity IV Birthing Bed

Product Description

- Overall Length: **90" (228.6 cm)**
- Length from Roller Bumpers to Break in Seat Section: **61 7/8" (157 cm)**
- Length from Break in the Seat Section to End of Bed: **28 1/8" (71 cm)**
- Maximum Width (siderails stored): **36" (91 cm)**
- Maximum Width (siderails up): **39" (99 cm)**
- Maximum Headboard Height: **52" (132 cm)**
- Maximum Siderail Height (without mattress): **14 3/4" (37.47 cm)**
- Minimum Under-Bed Clearance: **5" (13 cm)**
- Mattress Size: **35" x 78" (89 cm x 198 cm)**
- Caster Size: **6" (15.2 cm) or 8" (20 cm)**
- Maximum Bed Weight: **480 lb. (218 kg)**
- Head Section Inclination (maximum): **65°**
- Seat Section Inclination (maximum): **15°**
- Bed Height Range: **18" (46 cm) to 34" (86 cm)**
- Bed Height Range (with mattress): **22" (56 cm) to 38" (97 cm)**
- Trend-Like Position (maximum): **8°**
- Foot Section Weight Capacity (maximum): **400 lb. (181 kg)**
- Head Section Lift Capacity (maximum): **200 lb. (91 kg)**
- Maximum Patient Weight: **500 lb. (226.8 kg)**

