

Hill-Rom P881 GPS Stretcher

Product Description

- Overall Length: **84" (213.4 cm)**
- Height Overall –Low Position (without mattress): **23 ½" (59.7 cm)**
- Height Overall – High Position (without mattress): **38 ½" (97.8 cm)**
- Total Hi-low Travel: **15" (38.1 cm)**
- Overall Width – Siderail Stored: **28 ½" (72.4 cm)**
- Overall Width – Siderails Extended: **33 ½" (85.1 cm)**
- Siderail Length: **56 ¾" (144.1 cm)**
- Clearance Inside Siderails (width): **28" (71.1 cm)**
- Siderail Height Fully Up: **14 ½" (36.8 cm)**
- Siderail Height Intermediate: **12 ½" (31.7 cm)**
- Sleeping Surface: **26 ¼" x 78 ¼" (66.7 cm x 198.8 cm)**
- Trendelenburg Angle Inclination: **11°**
- Reverse Trendelenburg Inclination: **10°**
- Head Section Inclination – Assisted: **80°**
- Thigh Section Inclination: **40°**
- Maximum Stretcher Weight: **288 lb. (130.6 kg)**
- Maximum Patient Weight: **500 lb. (226.8 kg)**

***This model has a pneumatically operated head section, aluminum head pan cover, and non-adjustable foot and thigh section.**

