

# Stryker 1231 Renaissance Stretcher

## Product Description

- Overall Length: **84" (213.4 cm)**
- Overall Width (siderails up): **31 1/2" (80 cm)**
- Overall Width (siderails down): **30" (76.2 cm)**
- Height Overall / High Position: **35 1/2" (90.2 cm)**
- Height Overall / Low Position: **22" (55.9 cm)**
- Patient Surface: **26" x 76" (66 cm x 193 cm)**
- Siderails: **14" x 54" (35.6 cm x 137.2 cm)**
- Caster Size: **8" (20.3 cm)**
- Knee Gatch: **0° – 30°**
- Backrest: **0° – 90°**
- Trendelenburg/Reverse Trendelenburg: **±18"**
- Maximum Patient Weight: **500 lb. (226.8 kg)**

