

Stryker 2030 Epic II ICU Bed

Product Description

- Overall Length: **91" (231.1 cm)**
- Overall Width (siderails up): **42 ½" (108 cm)**
- Overall Width (siderails down): **36" (91.4 cm)**
- Height Overall – High Position: **32 ½" (82.5 cm)**
- Height Overall – Low Position: **18" (45.7 cm)**
- Caster Size: **6" (15.2 cm)**
- Knee Gatch Angle: **0° – 30°**
- Backrest: **0° – 90°**
- Trendelenburg/Reverse Trendelenburg: **±14°**
- Patient Surface: **35" x 84" (88.9 cm x 213.4 cm)**
- Patient Weight: **500 lb (226.8 kg)**

