

Stryker 3002 Secure II with Zoom

Product Description

- Overall Length: **93" (236.2 cm)**
- Overall Width (siderails up): **42 ½" (107.9 cm)**
- Overall Width (siderails down): **40" (102 cm)**
- Height Overall – High Position: **30" (76.2 cm)**
- Height Overall – Low Position: **16" (30.6 cm)**
- Caster Size: **6" (15.2 cm)**
- Fowler: **0° – 60°**
- Knee Gatch: **0° – 40°**
- Trendelenburg / Reverse Trendelenburg: **±12°**
- Retraction: **14" (35.6 cm)**
- Patient Surface: **35" x 84" (88.9 cm x 213.4 cm)**
- Patient Weight: **500 lb. (226.8 kg)**

